

## VERICHEM PRODUCTS FOR YOUR ANALYZER (2020)

Abbott Diagnostics:	ARCHITECT (all models)	All multi-level products compatible, Level (F) optional
	ALINITY (all models)	All multi-level products compatible, Level (F) optional
Beckman Coulter:	AU SYSTEMS (all models)	All multi-level products compatible, Level (F) optional
	DxC 700 AU	All multi-level products compatible, Level (F) optional
	UNICEL DxC (all models)	All multi-level products compatible, Level (F) optional
Roche Diagnostics:	AVL 9180 / 9181	9200 Electrolyte Standard Kit recommended
	COBAS C (all models)	All multi-level products compatible, Level (F) optional
	COBAS INTEGRA (all models)	All multi-level products compatible, Level (F) optional
Siemens Healthcare:	DIMENSION (RXL/EXL/Xpand )	All multi-level products compatible, Level (F) optional
	ATELLICA / ADVIA (all models)	All multi-level products compatible, Level (F) optional
	VISTA (all models)	All multi-level products compatible, Level (F) optional 9200 Electrolyte Standard Kit required for VISTA only
Advanced Instruments:	BR2 BILIRUBIN STAT	9450 Bilirubin, 9436 Bilirubin, Level F and 9459 BR2 Bilirubin 9456 Tru-Zero Bilirubin
Alfa Wassermann:	ACE Alera	All multi-level products compatible, Level (F) optional
EKF Diagnostics (Stanbio):	Altair / Excel	All multi-level products compatible, Level (F) optional
Horiba Medical:	PENTRA (all models)	All multi-level products compatible, Level (F) optional
		9240 Serum ISE Standard Kit Required
Instrumentation Labs:	GEM PREMIER (all models)	9500 MP Chemistry Reference and 9200 Electrolyte Kits
		9670 Ethanol and 9210 Carbon Dioxide Standard Kits
	GEM OPL & PCL PLUS	9500 MP Chemistry Reference and 9200 Electrolyte Kits
		9670 Ethanol and 9210 Carbon Dioxide Standard Kits
Ortho Clinical:	VITROS	9200 Electrolyte Standard Kit only
Medica Corporation:	EasyLyte	9240 Serum ISE Standard Kit only
	EasyRA	All multi-level products compatible, Level (F) optional
		9240 Serum ISE Standard Kit Required