

VERICHEM PRODUCTS FOR YOUR ANALYZER (2019)

Abbott Diagnostics:	ARCHITECT (all models)	All multi-level products compatible, Level (F) optional
	ALINITY (all models)	All multi-level products compatible, Level (F) optional
Beckman Coulter:	AU SYSTEMS (all models)	All multi-level products compatible, Level (F) optional
	DxC 700 AU	All multi-level products compatible, Level (F) optional
	UNiCEL DxC (all models)	All multi-level products compatible, Level (F) optional
Roche Diagnostics:	AVL 9180 / 9181	9200 Electrolyte Standard Kit only
	COBAS C (all models)	All multi-level products compatible, Level (F) optional
	COBAS INTEGRA (all models)	All multi-level products compatible, Level (F) optional
Siemens Healthcare:	DIMENSION (all models)	All multi-level products compatible, Level (F) optional
	ATELLICA (all models)	All multi-level products compatible, Level (F) optional
	ADVIA (all models)	All multi-level products compatible, Level (F) optional
Advanced Instruments:	BR2 BILIRUBIN STAT	9450 Bilirubin, 9436 Bilirubin, Level F and 9459 BR2 Bilirubin 9456 Tru-Zero Bilirubin
Alfa Wassermann:	ACE Alera	All multi-level products compatible, Level (F) optional
EKF Diagnostics (Stanbio):	Altair / Excel	All multi-level products compatible, Level (F) optional
Horiba Medical:	PENTRA (all models)	All multi-level products compatible, Level (F) optional
Instrumentation Labs:	GEM PREMIER (all models)	9500 MP Chemistry Reference and 9200 Electrolyte Kits 9670 Ethanol and 9210 Carbon Dioxide Standard Kits
	GEM OPL & PCL PLUS	9500 MP Chemistry Reference and 9200 Electrolyte Kits 9670 Ethanol and 9210 Carbon Dioxide Standard Kits
Ortho Clinical:	VITROS	9200 Electrolyte Standard Kit only
Medica Corporation:	EasyLyte	9200 Electrolyte Standard Kit only
	EasyRA	All multi-level products compatible, 9200 Electrolyte Required Level (F) optional
Randox Laboratories:	RX Series (all models)	All multi-level products compatible, Level (F) optional