

## RECOMMENDED VERICHEM PRODUCTS FOR YOUR ANALYZER

Beckman Coulter:	AU Systems (all models)	All products compatible - 9240 required for ISE
	DxC / LX / CX (all models)	All products compatible - 9240 required for ISE
Roche:	AVL 9180 / 9181	9200 Electrolyte Standard Kit only
	COBAS c501, c311, c411	All products compatible - 9240 required for ISE
	COBAS INTEGRA (all models) with direct or indirect ISE	All products compatible - 9240 required for ISE
	COBAS MIRA	All products compatible - 9240 required for ISE
	HITACHI (all models)	All products compatible - 9240 required for ISE
Siemens:	DIMENSION (all models)	All products compatible - 9240 required for ISE
	ADVIA	All products compatible - 9240 required for ISE
Abbott:	ARCHITECT	All products compatible - 9240 required for ISE
	AEROSSET	All products compatible - 9240 required for ISE
Elan:	ATAC	All products compatible - 9240 required for ISE
Medica:	EasyLyte	9200 Electrolyte Standard Kit only
	EasyRA	9500 MP Chemistry Reference and 9200 Electrolyte Kits only
Randox:	DAYTONA	All products compatible - 9240 required for ISE
Johnson & Johnson:	VITROS	9200 Electrolyte Standard Kit only
	FUSION	9200 Electrolyte Standard Kit only
Alfa Wassermann:	ACE	All products compatible - 9240 required for ISE
Advanced Instruments:	BR2 BILIRUBIN STAT	9450 Bilirubin, 9436 Bilirubin, Level F and 9459 BR2 Bilirubin
Stanbio:	SIRRUS	All products compatible - 9200 or 9240 required for ISE
Instrumentation Labs:	GEM PREMIER (all models)	9500 MP Chemistry Reference and 9200 Electrolyte Kits 9670 Ethanol and 9210 Carbon Dioxide Standard Kits
	GEM OPL & PCL PLUS	9500 MP Chemistry Reference and 9200 Electrolyte Kits 9670 Ethanol and 9210 Carbon Dioxide Standard Kits